


nebula

ENGLISH • FRANÇAIS • DEUTSCH • ESPAÑOL • PORTUGUÊS


A

10


11

12


13


14


15

21

20

19

16


17

22









23

18

do not disassemble these two parts
 ne pas démonter ces 2 parties
 die beiden teile nicht auseinanderbauen
 no desmontar estas dos partes
 não desmontar estas duas partes

B**C****D****E****F****G****H****J****I**

 Please read the instructions and the important notes (1.3) carefully before using the device.

1 - IMPORTANT INDICATIONS

1.1 - APPLICATION

NEBULA is an aerosol therapy device for administering drugs in spray form.

NEBULA incorporates a series of innovative solutions that allow rapid and effective therapeutic treatment of the upper and lower airways. The NEBULA SPACER kit (nebuliser + spacer mouth mask) or the paediatric mask are effective in treating disorders of the lower airways (including bronchial asthma, chronic bronchitis, pulmonary emphysema and cystic fibrosis).

The RINOWASH nasal douche is used for treating infections of the upper airways and for physiological flushing (with physiological or hypertonic solution or spa water) and/or therapeutic flushing (with drugs) of the nasal cavities.

 The administration of drugs requires a medical prescription specifying the type of drug, the dosage and the length of treatment.

This aerosol therapy device and its accessories have been designed and made in Italy in accordance with EU Directive 93/42/EC regarding medical devices (and subsequent updates) and bear the CE mark.

This product complies with the restriction on the use of certain hazardous substances in electric and electronic devices (European Directive RoHS 2011/65/EU).

Device suitable for continuous operation.

1.2 - SYMBOLS USED

Class II device  Device with type BF applied part  Important - read the instructions carefully 

Alternating current  Device ON "I"  Device OFF "O" 

The device must be dismantled at the end of its lifecycle and disposed of in special bins. Do not dispose of in ordinary town waste bins (EU Directive 2002/96/EC) 

Device protected against finger access to hazardous parts but not protected against the ingress of liquids (IP20)

Keep the device dry  Manufacturer 

Does not contain phthalates  No traces of natural latex 

1.3 - IMPORTANT NOTES

- Read the instructions carefully and keep them in a safe place for future reference.
- Do not use the device or accessories for anything other than the intended purpose.
- Children must only be allowed to use the device under the supervision of an adult who has read and understands the instructions.
- Some components are small enough to swallow and may lead to suffocation.
- The power cable and connecting pipe are long enough to pose a strangulation hazard.
- Read the list of components in the table in section 7 to see if there are any components that have caused an allergic reaction in the past.
- Do not use essential oils (e.g. with menthol or eucalyptus) as they are not compatible with the material of which the nebuliser is made.
- During the session, the device must stand on a firm flat surface.
- Never cover the air vents during operation.
- Keep well away from sources of heat.
- Do not touch the plug or the electrical components wet or with damp hands.
- Device protected against finger access to hazardous parts but not protected against the ingress of liquids (IP20).
- Do not pour water or other liquids onto the device. Prevent water and other liquids from getting in through the air vents.
- Do not use in the bath or shower.
- If the device should fall into liquid, unplug before removing it. Do not use the device. Send it to an approved service centre for an overhaul.
- Beware of compressed air coming out of the air outlet as this may cause a hazard; it must only be used to supply the nebuliser.
- The device is not suitable for use in the presence of anaesthetic mixtures with air, oxygen or nitrogen protoxide.
- Performance of the device is only guaranteed if all the accessories used are original ones. Only original spare parts must be used.
- The manufacturer declares that the device has a useful life of at least 5 years.

1.4 - DESCRIPTION OF THE DEVICE

The NEBULA aerosol therapy device is made up of a compressor driven by an electric motor in the basic unit, a NEBULA SPACER kit (nebuliser + spacer mouth mask), a paediatric mask (which can be mounted on the nebuliser as an alternative to the spacer mouth mask), a RINOWASH nasal douche and a flexible tube connecting the nebuliser/nasal douche to the basic unit.

1.5 - PATIENT GROUPS

The NEBULA aerosol therapy device provides highly effective inhalation therapy for patients of any age, from young children to the elderly.

1.6 - CONTENTS

The package contains (Fig. A):

the basic unit (1) with the power cable attached, a cable storage compartment (2), an accessory storage compartment (3), an ON/OFF switch (4), a filter slot (5), a fuse holder (6), an air intake (7), and a cover (8) with a support (9) for the nebuliser.

It also contains:

- a Nebula Spacer kit containing the following:
 - a nebuliser (upper section (10), atomiser (11), lower section (12))
 - a spacer mouth mask (mask body (13), fitting with discharge valve (14))
- a paediatric mask (15)
- a Rinowash nasal douche made up of the following:
 - a base with a button (16), a connecting element for the tube (17), a base cap (18) and a central tank (19)
 - atomiser (20)
 - an outer bell (21)
- a connecting tube (22)
- spare air filters (23)

2 - PREPARATION

2.1 - PREPARING THE BASIC UNIT

Before use, make sure the voltage indicated on the rating plate underneath the device corresponds to the mains voltage and the wiring system complies with the current regulations.
Do not use if the cable is damaged.
Check that the casing is integral and there is no sign of dents or damage.

Make sure the switch (4) is in position 0 (OFF). Insert the electric plug in the mains socket, press position I (ON) on the switch and check that air comes out of the air intake (7). Then press position 0 (OFF) on the switch.

2.2 - PREPARING THE NEBULA SPACER KIT

Please read section 3.2 before using the device.

Proceed as follows to use the Nebula Spacer:

- Make sure the atomiser (11) is correctly inserted in the central sleeve in the lower section (12) of the nebuliser
- The lower section (12) has a graduated scale that shows the approximate amount of liquid inside. The purpose is merely to give an approximate indication of the contents. Use a graduated syringe to get an accurate dose. Pour the liquid into the tank (Fig. B), keeping to the dosage specified by your doctor.

Readings must be taken on the graduated scale with the atomiser mounted on the lower section.

- Screw the upper section (10) of the nebuliser onto the base.
- Connect the fitting (14) of the spacer mouth mask to the upper section (10) of the nebuliser (Fig. C).
- Connect the lower end of the tube (22) to the nebuliser and the other end to the air intake (7).
- Press position I (ON) on the switch (4).
- During treatment, hold the nebuliser upright in your hand to enable it to work efficiently.
- If the above operations have been performed correctly, atomised spray will come out of the mask.

The spacer mouth mask (13) must be placed gently over the mouth – the top of the mask rests underneath the nose. For the treatment to be effective, breathe in through the mouth and out through the nose. If the user breathes out through the mouth, the valve on the fitting (14) releases the exhaled air (Fig. D).

Do not obstruct the valve on the fitting during use (Fig. D).

When you have used up all the liquid or wish to interrupt treatment, press (0) (OFF) on the switch and place the nebuliser back on its support (9).

A small amount of liquid may remain in the nebuliser at the end of the session but this has no importance.

If a small amount of moisture remains inside the air tube at the end of the session, leave the device switched on and allow it to dry by ventilation.

The presence of liquid inside the tube is the result of water vapour forming under certain environmental conditions.

2.3 - USING THE PAEDIATRIC MASK

Please read section 3.3 before using the device.

To use the paediatric mask, fit it onto the upper section of the nebuliser, in the same way as for the spacer mask (Fig. C).

With children, the paediatric mask must be placed gently over the mouth and nose. The treatment consists of inhaling the atomised liquid slowly and deeply. The exhaled air is released through vents in the sides of the mask.

During use, do not obstruct the air vents at the sides of the mask.

2.4 - USING THE RINOWASH NASAL DOUCH

Please read section 3.4 before using the device.

Proceed as follows to use the Rinowash:

Pour up to 15 ml of liquid into the tank (Fig. E), reading the amount on the graduated scale.

Readings must be taken on the graduated scale only when the atomiser is not mounted on the tank.

- Fit the atomiser (20) onto the sleeve of the tank. Make sure it clicks into position.
- Fit on the outer bell (21), turning slightly, to cover the seal.
- Attach one end of the connecting tube (22) to the Rinowash connector (17) and the other end to the air intake (7).
- Press I (ON) on the switch (4) to start up.
- Place the end of the Rinowash in one nostril (Fig. F).
- Press the button on the base (16) to generate a spray.
- Therapeutic spray is delivered through the hole in the outer bell (21), and residual fluid collects in the chamber between the outer bell and the tank (Fig. G).
- During treatment, breathe normally through the nose. If mucus forms in the nostril, blow it out. The space between the outer bell and the tank acts as a chamber for collecting the fluid that drains out of the nose.
- Treat the nostrils alternately, as required.

Atomisation is very rapid and the session only lasts a few minutes.

Make sure all the parts are connected firmly together. Check that the ends of the tubes are inserted properly and the components of the Rinowash are mounted correctly.

If the device is not assembled correctly, atomisation of the liquid may be affected. If this happens, rotate the atomiser (20).

As stated under Important Notes, some liquids are not compatible with the materials of which the nebuliser, Rinowash and spacer mouth mask are made. They may interact and cause damage (the materials are shown in section 7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS).

3 - CLEANING, DISINFECTION AND STERILISATION

Before proceeding, unplug the device and disconnect the nebuliser, Rinowash and tube.

If it is necessary to clean, disinfect or sterilise this medical device using a method other than that given in this section, contact the manufacturer first.

3.1 - CLEANING, DISINFECTING AND STERILISING THE BASIC UNIT

3.1.1 - CLEANING THE BASIC UNIT

The basic unit (1) must not be cleaned in water as it could get inside and damage the electrical components, which would cause a safety hazard. Clean the basic unit (1) with a soft damp cloth, with a small amount of mild (non-abrasive) detergent, if necessary.

3.1.2 - DISINFECTING THE BASIC UNIT

The basic unit can be disinfected using a soft sponge or cloth and some cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions.

Make sure the solution is compatible with the material of which the basic unit is made. Refer to section 7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS.

3.1.3 - STERILISING THE BASIC UNIT

The basic unit (1) must not be sterilised.

3.2 - CLEANING, DISINFECTING AND STERILISING THE NEBULA SPACER KIT

To prevent the risk of microbial contamination, the nebuliser and the spacer mouth mask must be cleaned and disinfected, and if necessary sterilised, after each treatment and before use, following the instructions carefully.

Before cleaning, disinfecting or sterilising any of the components, disassemble them as follows:

- Remove the tube (22).
- Remove the spacer mouth mask (13, 14) from the nebuliser.
- Detach the fitting (14) from the mask (13).
- Unscrew the upper section (10).
- Remove the atomiser (11).

3.2.1 - CLEANING THE NEBULA SPACER KIT

The disassembled nebuliser (components 10, 11 and 12) and the disassembled spacer mouth mask (13, 14) can be washed in water at an approximate temperature of 40°C with the addition of detergent (refer to the manufacturer's instructions), then rinsed in running water. When all the parts have been cleaned, dry them with a soft cloth and place in the storage compartment.

3.2.2 - DISINFECTING THE NEBULA SPACER KIT

After washing, all the components of the nebuliser (10, 11 and 12) and the spacer mouth mask (13, 14) can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions. The components of the nebuliser and the spacer mask can be disinfected by boiling in a saucepan of water for a maximum of 10 minutes. Make sure the components are not in direct contact with the bottom of the saucepan.

Boiling may alter the colour or appearance of the components, but this will not affect operation of the device.

3.2.3 - STERILISING THE NEBULA SPACER KIT

All the components of the nebuliser (10, 11 and 12) and the spacer mouth mask (13, 14) can be sterilised in an autoclave at 121°C for at least 15 minutes (but no more than 30 minutes), for up to 20 times.

Sterilisation is recommended before each session when the device is used in a hospital environment or by more than one patient.

3.3 - CLEANING, DISINFECTING AND STERILISING THE PAEDIATRIC MASK

To prevent the risk of microbial contamination, the paediatric mask (15) can be cleaned and disinfected, and if necessary sterilised, after each treatment and before use, following the instructions carefully.

3.3.1 - CLEANING THE PAEDIATRIC MASK

The mask can be washed in water at an approximate temperature of 40°C with the addition of detergent (refer to the manufacturer's instructions), then rinsed in running water. When the mask has been cleaned, dry it with a soft cloth and place in the storage compartment.

3.3.2 - DISINFECTING THE PAEDIATRIC MASK

After washing, the mask can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite), following the manufacturer's instructions.

The paediatric mask (15) must not be boiled

3.3.3 - STERILISING THE PAEDIATRIC MASK

The paediatric mask (15) must not be sterilised hot or in an autoclave.

3.4 - CLEANING, DISINFECTING AND STERILISING THE RINOWASH NASAL DOUCHE

To prevent the risk of microbial contamination, the Rinowash can be cleaned and disinfected, and if necessary sterilised, after each treatment and before use, following the instructions carefully.

Before cleaning, disinfecting or sterilising any of the components, disassemble them as follows:

- Remove the tube (22).
- Remove the outer bell (21) by turning gently.
- Remove the atomiser (20) by pressing on the tab (Fig. H).
- Remove the base cap (18).

Do not use a microwave oven for cleaning, disinfecting or sterilising.

The seal at the base of the tank (19) must not be removed. The main section and the tank (19) form a single component and cannot be disassembled.

3.4.1 - CLEANING THE RINOWASH

The Rinowash can be washed in hot water using non-abrasive detergent and rinsed in running water. When all the parts have been cleaned, dry them with a soft cloth, reassemble them and place in the storage compartment.

3.4.2 - DISINFECTING THE RINOWASH

After washing, all the components of the Rinowash can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions. N.B. Make sure the disinfectant solution is compatible with the materials used. The components can be disinfected by boiling in a saucepan of water for a maximum of 10 minutes. Make sure the components are not in direct contact with the bottom of the saucepan.

Boiling may alter the colour or appearance of the components, but this will not affect operation of the device.

3.4.3 - STERILISING THE RINOWASH

All the components of the Rinowash can be sterilised in an autoclave at 121°C for at least 15 minutes (but no more than 30 minutes), for up to 20 times.

3.5 - CLEANING, DISINFECTING AND STERILISING THE TUBE

The tube (22) needs to be cleaned after each session. Detach it from the nebuliser and leave the device on until all signs of moisture have been eliminated. If the tube is particularly dirty, it needs to be replaced.

 The tube (22) must not be sterilised hot or in an autoclave or boiled.

4 - MAINTENANCE

The NEBULA aerosol therapy device does not require any particular maintenance. The compressor does not need to be lubricated as the piston is the dry-operating type.

A special external filter keeps out dust and impurities, and should be replaced periodically.

 Always unplug from the mains before carrying out the following operations.

4.1 - REPLACING THE FUSE (Fig. I)

If the device does not work even when plugged in, it may be necessary to change the fuse.

Unplug the device.

Use a small screwdriver (not included) to remove the fuse holder, which is located next to the switch.

Remove the fuse and replace it with a new one having the same features (see section 6 - TECHNICAL DATA).

Insert the fuse holder in the slot.

Plug in and check to see if it switches on.

4.2 - REPLACING THE AIR FILTER (Fig. J)

The air filter needs to be replaced if it gets dirty, and once a year if the device is used frequently.

To replace the air filter (5), remove the cap using a screwdriver and insert a new one.

The filter cannot be washed or regenerated. It must always be replaced with a new one.

4.3 - MAINTENANCE OF THE ACCESSORIES

When cleaning, disinfecting or sterilising the components, you can check them visually and determine when they need to be replaced. Repeated sterilisation may affect the materials of which the accessories are made, so the components should be inspected regularly. The accessories must be replaced when there are clear signs of deterioration.

5 - PROBLEMS, POSSIBLE CAUSES AND SOLUTIONS

TROUBLE-SHOOTING CHART

Problems	Possible causes	Solutions
The device does not work.	<ul style="list-style-type: none">- The plug is not inserted in the socket.- The switch (4) is in position 0 (OFF).- The fuse has blown.- The thermostat has cut in.	<ul style="list-style-type: none">- Insert the plug in the socket and press position I (ON) on the switch (4).- Press position I (ON) on the switch (4).- Replace the fuse (see 4.1).- Clear the air vents.
The nebuliser does not generate a spray.	<ul style="list-style-type: none">- The device is switched off.- There is no liquid in the nebuliser.- Too much liquid has been poured into the nebuliser.- The nebuliser has not been assembled correctly.- The air outlet in the base (12) is blocked.- The air tube (22) is not connected securely.	<ul style="list-style-type: none">- Switch on the device.- Pour liquid into the nebuliser.- Reduce the quantity of liquid in the nebuliser to below the maximum level.- Check that the atomiser (11) is properly inserted in the base of the nebuliser (12). Screw the upper section (10) onto the base and check that air comes out of the little hole at the top of the base.- Unblock the hole.- Connect the tube (22) securely to the air intake on the device (7) and the one on the nebuliser.
The nebulised liquid does not reach the mask or there is difficulty exhaling (problems due to the spacer mask).	<ul style="list-style-type: none">- The spacer mask (13, 14) has not been assembled correctly.- The valve (Fig. D) has jammed.	<ul style="list-style-type: none">- Assemble the spacer mask (13, 14) correctly.- Release the valve.

Problems	Possible causes	Solutions
The Rinowash does not deliver a spray.	<ul style="list-style-type: none"> - The device is switched off. - There is no liquid in the tank. - The air tube (22) is not connected securely or the Rinowash has not been assembled correctly. 	<ul style="list-style-type: none"> - Switch on the device. - Pour liquid into the tank. - Connect the tube (22) securely to the air intake on the device (7) and the one on the Rinowash (17). Rotate the atomiser (20) slightly. Make sure all the parts have been assembled correctly.
Liquid leaks from the base of the bell.	<ul style="list-style-type: none"> - The outer bell (21) has not been fitted on properly and the seal does not hold. 	<ul style="list-style-type: none"> - Fit on the outer bell (21) correctly, turning slightly. Make sure the whole seal is covered.

If you have failed to solve the problem after performing the above checks and the device still does not work, have it checked by an authorized service centre.


6 - TECHNICAL DATA

Supply voltage	See rating plate on the device
Operation	Continuous
Maximum air pressure	2.7 bar
Maximum air flow	12 l/min
Dynamic air flow with nebuliser	7 l/min
Maximum input	See rating plate on the device
Safety fuse	See rating plate on the device
Dimensions of the basic unit	230x180x185H mm
Weight of the basic unit	3 kg
Noise level at 1 m	57 dBA (to UNI EN 13544 -1)
Minimum nebuliser filling volume	1ml
Maximum nebuliser filling volume	8 ml
Minimum Rinowash filling volume	1,5 ml
Maximum Rinowash filling volume	15 ml
Aerosol output with nebuliser (*)	0.75 ml (**)
Minimum aerosol output rate with nebuliser (**)	0.30 ml/min (**)
MMAD with nebuliser, graph 1 (*) (***)	1.9 µm (**)
MMAD with Rinowash micronized nasal douche, graph 2 (***)	18 µm (**)


The MMAD values shown are measured using a MACH 2 API AEROSIZER.

The y-axis of the graphs represents the percentage of volume of atomised solution with sizes lower than the corresponding dimensional value expressed in µm on the x-axis.

GRAPH 1

GRAPH 2

(*) For a quantity of 2 ml.

(**) The values shown refer to the use of physiological solution (0,9% NaCl) and may vary with the drug used. The values given do not apply to drugs in suspension or high-viscosity drugs. The manufacturer of the drug can provide the necessary details.

(***) 50% by volume of the atomised particles is less than the MMAD (Median Mass Aerodynamic Diameter). This value is indicative of the nebuliser's ability to produce particles suitable for a specific therapeutic treatment.

7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS

At the end of their useful life, the device and its components must be disposed of in accordance with the current regulations and environmental laws. If no legal obligations exist, the various materials must be sorted and disposed of separately.

The table below shows the materials of which the various components are made.

None of the materials used contains phthalates or traces of natural Latex.

COMPONENT	MATERIAL
Basic unit (1)	Waste from Electric and Electronic Equipment (WEEE)
External casing of the basic unit	ABS/Polystyrene
Cover of the basic unit	Polystyrene
Nebuliser (10, 12)	Polycarbonate
Atomiser (11)	Nylon
Mask body (13) and fitting (14)	Silicone and polycarbonate
Paediatric mask (15)	Non-toxic PVC (phthalate-free polyvinyl chloride)
Rinowash body (16, 17, 19)	Polycarbonate (some metal parts)
Rinowash base cap (18)	Elastomer
Rinowash atomiser (20)	Polycarbonate
Rinowash outer bell (21)	Polycarbonate and elastomer
Connecting tube (22)	Non-toxic PVC (phthalate-free polyvinyl chloride)
Box	Cardboard
Instructions	Paper
Accessory bag	Low-density polyethylene (LDPE)

8 - ELECTROMAGNETIC COMPATIBILITY

This device complies with the requirements of CEI EN 60601-1-2:2010 on the electromagnetic compatibility of medical devices.

Compliance with electromagnetic compatibility standards does not guarantee total immunity of the product. Some devices (mobile phones, pagers, etc.) may interrupt operation if used near medical devices.

RECOMMENDED DISTANCES BETWEEN RADIO-COMMUNICATION DEVICES

The NEBULA aerosol therapy device is designed to operate in an electromagnetic environment where radiated RF disturbance is under control. The customer or user of the NEBULA aerosol therapy device can help prevent electromagnetic interference by leaving a minimum distance between mobile and portable RF communication devices (transmitters) and the NEBULA aerosol therapy device. The table below shows the recommended distances for the NEBULA aerosol therapy device based on the maximum output of the radio communication devices.

Nominal maximum transmitter output [W]	Separation distance [m] at the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \times \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \times \sqrt{P}$
0.01	0,12	0,12	0,23
0.1	0,38	0,38	0,73
1	1,2	1,2	2,33
10	3,8	3,8	7,3
100	12	12	23

The recommended separation distance d in metres [m] for transmitters with a nominal maximum transmitter output not shown in the table can be calculated using the equation applied to the transmitter frequency, where P is the nominal maximum transmitted output in Watt [W] declared by the manufacturer.

Note 1. For 80 MHz and 800 MHz, apply the distance for the higher frequency interval.

Note 2. These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and persons.

9 - OPERATING AND STORAGE CONDITIONS

9.1 - OPERATING CONDITIONS

- Temperature range +5°C to +40°C
- Relative humidity 15% to 95%
- Atmospheric pressure 700 hPa to 1060 hPa

9.2 - STORAGE CONDITIONS

- Temperature range -25°C to +50°C
- Relative humidity less than 93%
- Atmospheric pressure 700 hPa to 1060 hPa

10 - SPARE PARTS

Only original spare parts must be used.

11 - WARRANTY CONDITIONS

The manufacturer is liable for the safety, reliability and performance of the device, provided that it is used in accordance with the instructions and for the intended purpose only, and any repairs are carried out by the manufacturer or an approved service centre.

The manufacturer declines all liability and the warranty will be invalidated if the user fails to follow the instructions and the important notes contained herein.

12 - TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated accordingly to include any new or changed features.

If the booklet accompanying the device gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

13 - COPYRIGHT

The information contained in this booklet must not be used for any purpose other than the reason for which it is provided.

This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

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