MANUEL D'UTILISATION • USER GUIDE • GEBRAUCHSANLEITUNG • MANUAL DE USO • ISTRUZIONI D'USO • INSTRUÇÕES DE USO • BRUKSANVISNING
GEBRUIKSAANWIJZING • INSTRUKCJA OBSŁUGI • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ • ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ • UŽIVATELSKÝ MANUÁL • 使用说明 • 使用説明書


Respireo Primo F • non vented

MASQUE BUCCO-NASAL • FACIAL MASK • MUND-NASEN-MASKE • MASCARILLA ORO-NASAL • MASCHERA ORO-NASALE
MÁSCARA ORO-NASAL • HELMASK • MOND-NEUSMASKER • MASKA USTNO-NOSOWA • ΠΙΝΙΚΗ ΣΤΟΜΑΤΙΚΗ ΜΑΣΚΑ
ПОТО-НОСОВАЯ МАСКА • NOSO-ÚSTNÍ MASKA • 口鼻面罩 • 鼻口腔マスク


| | |
|-----------------|-------|
| FRANÇAIS..... | p. 7 |
| ENGLISH..... | p. 10 |
| DEUTSCH..... | p. 13 |
| ESPAÑOL..... | p. 17 |
| ITALIANO..... | p. 20 |
| PORTUGUÊS..... | p. 24 |
| SVENSKA..... | p. 27 |
| NEDERLANDS..... | p. 30 |
| POLSKI..... | p. 34 |
| ΕΛΛΗΝΙΚΑ..... | p. 37 |
| РУССКИЙ..... | p. 41 |
| ČEŠTINA..... | p. 45 |
| 中文..... | p. 48 |
| 日本語..... | p. 51 |

A


B


1


2





3



4



B

C

D

1 - UTILISATION

Respireo Primo F Non Vented est un masque bucco-nasal sans fuites calibrées (Non Vented) à patient unique.

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1 - APPLICATION

Respireo Primo F Non Vented is a non-vented facial mask for use by a single patient.

Respireo Primo F Non Vented is designed for use by patients with controlled non-invasive positive pressure ventilation (NPPV).

Respireo Primo F Non Vented is suitable for adult patients (>30Kg).

⚠ 2 - IMPORTANT NOTES

- 2.1 The mask must only be used with a ventilator equipped with an active exhalation valve or a respiratory circuit with an exhalation valve.
- 2.2 The operating pressure range is 4-40 cmH₂O.
- 2.3 If you experience discomfort or irritation, or suffer an allergic reaction to any of the components of the mask, please report this to your doctor or respiratory therapist.
- 2.4 If oxygen is used with the ventilator, the supply must be turned off when the appliance is not in operation.
- 2.5 Do not smoke or light naked flames when using oxygen.
- 2.6 Do not use the mask in the presence of nausea or vomiting.
- 2.7 The mask must be cleaned and/or disinfected. For further details refer to section 7 - CLEANING AND DISINFECTING THE MASK.
- 2.8 Check the mask before use. If there is any visible damage, report it to the dealer.
- 2.9 The mask must only be used in accordance with these instructions. No other instructions are applicable.
- 2.10 Some of the components of the mask are small enough to be swallowed by children and must never be left unattended. An adult must always be present.

3 - COMPONENTS OF THE MASK (Fig. A)

a) mouth-and-nose cushion; b) main body; c) forehead support; d) forehead rest; e) adjusting knob; f) adjusting screw; g) cap; h) elbow; i) headgear; l) clip.

4 - HOW TO PUT ON THE MASK (Fig. B)

The mask comes fully assembled ready for use. Proceed as follows to put on the mask.

Turn the adjusting knob (e) until the forehead support (c) is fully extended (Fig. 1).

Hold the mask firmly in one hand and extend the headgear (i) with the other enough to pull it over your head (Fig. 2).

Place the mask gently on your face so that the cushion (a) fits snugly round the mouth and nose, and position the headgear (i) on your head (Fig. 3).

Clip the two clips (l) onto the main body (b) (Fig. 4).

Tighten the upper straps first, then the lower straps, until the mask fits on your face comfortably, without over tightening (Fig. 5).

Connect the mask to the ventilator circuit (Fig. 6)

After lying down, turn the adjusting knob (e) until the mask feels comfortable but still provides a good seal (Fig. 7).

Readjust the lower straps if necessary to eliminate any leaks, but make sure the mask is not too tight.

5 - HOW TO REMOVE THE MASK (Fig. C)

To remove the mask, press the button on the front (Fig. C) to release one of the clips (l). Lift the headgear over your head and remove the mask.

The other straps do not need to be loosened or disconnected.

6 - HOW TO DISASSEMBLE THE MASK

Before disassembling the mask, disconnect the ventilator circuit.

Make a note of the position of the adjusting knob (e) for reference next time.

Unclip the two clips (l) from the main body (b).

Undo the four straps and detach the headgear (i) from the main body (b) and the clips (l).

Pull the elbow (h) off the main body (b).

Remove the cushion (a) from the main body (b). Remove the cap

(g) from the main body (b).

Detach the forehead rest (d) from the forehead support (c).

7 - CLEANING AND DISINFECTING THE MASK

The mask must be cleaned and/or disinfected each time before use.

N.B. The headgear should be cleaned once a week.

CLEANING THE MASK

First wash your hands thoroughly.

Disassemble all the components (see 6 - HOW TO DISASSEMBLE THE MASK) and wash carefully in water at about 30°C using mild soap. Rinse thoroughly with water and leave in the air to dry.

The mask can be cleaned in a washing machine at 30°C (programme for delicate fabrics) but must not be spin- or tumble-dried. The mask can be washed in this way up to 20 times.

DISINFECTING THE MASK

The mask must be taken apart for disinfection (see 6 - HOW TO DISASSEMBLE THE MASK). Immerse the components in a cold solution of disinfectant (e.g. a solution containing up to 2% sodium hypochlorite), following the manufacturer's instructions.

IMPORTANT. Check that the disinfectant solution is compatible with the materials of which the mask is made (see 12 - DISPOSAL INSTRUCTIONS).

N.B. The headgear (i) must not be disinfected.

8 - PRECAUTIONS

Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents or aromatic oils. Such solutions would damage the product and reduce its life. Prolonged exposure to direct sunlight or heat may cause damage to the mask.

If any visible damage is noted (cracks, tears, etc.), the damaged part must be replaced.

Check carefully to make sure there are no soap deposits or foreign bodies inside the mask.

9 - HOW TO REASSEMBLE THE MASK

Insert the elbow (h) in the main body (b) of the mask, pressing gently.

Remount the forehead rest (d) on the forehead support (c), keeping the curved edge facing upwards.

Fit the grooved side of the cushion (a) into the frame of the main body (b).

Fit on the top first, then the sides of the mask, pressing the cushion in until the whole edge is up against the six stops on the main body (Fig. D).

Fit on the cap (g).

Pass the two lower straps of the headgear (i) through the slots in the clips (l).

Attach the two upper straps of the headgear (i) to the main body (b). Clip the two clips (l) back onto the main body (b).

Turn the adjusting knob (e) to the reference position you noted previously.

10 - PRESSURE/OXYGEN PORT

If additional oxygen or a pressure reading is required, use the port in the main body (b).

Open the cap (g) by lifting upwards, without removing it.

Push the pressure/oxygen tube firmly into the port.

11 - TROUBLE-SHOOTING CHART

| PROBLEMS | POSSIBLE CAUSES | SOLUTIONS |
|--|---|--|
| • No air or insufficient air comes out of the mask | <ul style="list-style-type: none"> • Flow generator off • Flow generator regulated wrongly • Patient circuit tube not inserted | <ul style="list-style-type: none"> • Switch on the flow generator • Regulate the flow generator • Connect the tube to the flow generator and mask |

| | | |
|---|--|-------------------------------|
| • The mask is noisy | • Cushion mounted incorrectly and air leak | • Mount the cushion correctly |
| • The mask does not fit the face snugly | • Mask applied incorrectly | • Reposition the headgear |

12 - DISPOSAL INSTRUCTIONS

The components of the mask are made of the following materials and must be disposed of in accordance with the applicable laws.

| COMPONENT | MATERIAL |
|------------------|--------------------|
| Main body | Polycarbonate |
| Forehead support | Polycarbonate |
| Adjusting knob | Polypropylene |
| Adjusting screw | Polycarbonate |
| Elbow | Polycarbonate |
| Forehead rest | Silicone |
| Cushion | Silicone |
| Cap | Silicone |
| Headgear | Nylon/polyurethane |
| Clip | Nylon |

13 - TECHNICAL DATA

| | |
|--|--|
| Normative references | CE marking in accordance with European Directive 93/42/EEC, class IIa and subsequent amendments EN ISO 14971, ISO 5356-1, ISO 10993-1, EN ISO 17510-2 |
| Operating pressure | 4-40 cmH ₂ O |
| Temperature range | 5°C - 40°C |
| Dead space volume inside the mask | 215 ml (Large) |

14 - LIMITED WARRANTY

Air Liquide Medical Systems S.r.l. warrants that your product shall be free from defects in material and workmanship for a period of 120 days from the date of purchaser.

If the product fails under conditions of normal use, Air Liquide Medical Systems S.r.l. will repair or replace, at its option the defective product or any of its components.

The manufacturer is liable for the safety, reliability and performance of the apparatus, provided that it is used in accordance with the instructions and for the intended purpose only, and that any repairs are made by the manufacturer or an approved service centre.

Throughout the warranty period keep the original packing.

Warranty claims on effective product must be made by the consumer contacting the dealer.

For further information on your warranty rights contact your dealer.

COPYRIGHT

The information contained in this booklet may not be used for anything other than the intended purpose. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated to ensure that they continue to be consistent with the characteristics of the devices marketed.

If the booklet accompanying the apparatus is mislaid, a replacement copy can be obtained from the Manufacturer by quoting the data shown on the label.

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AIR LIQUIDE MEDICAL SySTEMS S.p.A.