


# ***Monnal T75***

Tailor-made breathing



[www.airliquidemedicalsystems.com](http://www.airliquidemedicalsystems.com)

# Top-class performance

A ventilator operates in varied clinical situations specific to every patient. Monnal T75 provides a new solution due to its high level of performances and polyvalence. It combines simplicity of use and servicing and an excellent integration in all different hospital's environments.



**Optimal Breathing System**

**M**onnal T75 is one of a new generation of ventilators intended for hospital, to be used for critical, intensive and intermediate cares. It allows both invasive and non-invasive ventilation for patients from children to adults.

**Monnal T75** is the result of over 30 years' experience in our mechanical ventilation laboratories and of our collaboration with medical teams. It provides efficient management of the most acute pathologies.

**Monnal T75** contributes to increase the comfort of each patient by:

- leakage management and compensation by means of a cycle-to-cycle leak index,
- flexible trigger reactivity,
- excellent pressurization for providing significant reductions in efforts to be made by patients.



**Monnal High Quality**

**O**ur commitment is to provide you with efficient ventilation throughout **Monnal T75's** service life. The high quality of this ventilator relies more especially on:

- well-tested and proven components,
- materials selected for their resistance and their high thermal and acoustic capacity,
- quality requirements in its manufacture and finishing.

## Designed for handling the most acute pathologies

label  
observeurdudesign<sup>10</sup>

Monnal T75's design has been thought out perfectly for harmonious integration into hospital architecture and the unit's medical environment. It has also been designed for the serenity of nursing staff, patients and their families:

- relaxing colors,
- rounded forms,
- reassuring materials chosen with care.

The Observeur du design' Label was awarded to Monnal T75 by a judging panel from APCI (Agency for the Promotion of Industrial Creation). Awarded annually for more than 10 years, the Observeur du design Label rewards the most outstanding products resulting from collaboration between manufacturers and designers.

# monnal T75

## Complete patient monitoring

- Ventilation curves and loops choice
- 8 patient's parameters permanently measured and displayed
- Direct access to alarms settings

## Advanced functions

- Alarms history
- Ventilatory mechanic exploration
- Integrated nebulization



# Unique simplicity of use

Monnal T75 is a new but simple device. Its rapid, comfortable and reliable operation frees you completely for carrying out your health-care mission with patients.



*Universal Vision*



*Natural Interface*

**W**ith **Monnal T75** you benefit from a wide view for piloting your ventilator in all simplicity and calm thanks to:

- a tactile colors screen (10"4) for perfect visibility,
- a single page concept with direct access to:
  - set points,
  - alarms,
  - monitoring,
  - ventilation curves and loops,
- a highly structured page with automated operating and information-reading systems.

This system saves time and frees you for attending your patients.

**T**his single structured page concept, enables you to use a new generation ventilator with a single intuitive principle of operation:

- on-screen tactile selection,
- settings and validation made by code wheel.



## Designed for optimizing the time dedicated to each patient

# Integration adapted to the environment

Monnal T75 has been designed to cover all requirements related to patients' and nursing staff's environment: autonomy, compactness, mobility. It also meets new constraints in utilization costs, hygiene and sterilization.

**M**onnal T75 accompanies patients during intra-hospital transports thanks to:

- its autonomy in gas and electricity,
- the ergonomic design of its support unit.




*Easy Valve*

## Costs control

- An excellent capability / price ratio,
- simplified preventive maintenance that consists of a single annual performance check,
- medical team training optimized thanks to the extremely simple interface,
- compatible with the existing equipment pool, thanks to Monnal T75's polyvalence.

**D**aily servicing of **Monnal T75** is easy, rapid and reliable. Its innovative expiratory unit contains:

- only three components for even more ease in handling and sterilization,
- a single use or autoclavable use model,
- an ergonomic ejector button,
- simple click-locking.



## Designed to suit all medical environments



## Contact

Air Liquide Medical Systems  
*Parc de Haute Technologie*  
*6, rue Georges Besse*  
*92182 Antony Cedex, France*  
*Tel: +33 (0)1 40 96 66 00*  
*Fax: +33 (0)1 40 96 67 00*



Air Liquide's Healthcare business supplies medical gases & equipment to hospitals, medical services to patients at home and is active in the fight against nosocomial infections. Its 7,800 employees serve 5,000 hospitals and 500,000 patients, on five continents. Air Liquide is the world leader in gases for industry, health and the environment.