

monnal T30

Excellence serving well-being



Monnal_{T30} belongs to the last generation of Bi-Level pressure ventilators, created for patients who need the highest level of respiratory comfort.



Optimal Breathing System

- ⌚ 4 inspiratory triggers in order to naturally **fit** to any patient need
- ⌚ A recognized sensitivity of the triggers, in order to offer **unrivalled softness**
- ⌚ A **continuous** measure of leaks and tidal volume, regardless the mask or the pressure level be selected



Natural Interface

- ⌚ **Intuitive** to encourage a daily use
- ⌚ **Easy** to access to main functions



Homecare Adaptive Device

- ⌚ A **modern design** to discretely fit in patient's home



Monnal High Quality

- ⌚ A selection of components guided by an objective of **longevity** superior to the market average

Technical Specifications

Modes: CPAP, S, ST, ACMV, T
Inspiratory pressure: 4 to 20hPa (MonnalT20),
 4 to 30 hPa (MonnalT30)
Tins: 0.3 to 5s
Rise time: 200 to 600 ms/ 10cm H2O
Alarms: audible & visual
Sound level: 30dB
Power supply: 110 – 230VAC, 50-60Hz
Size: (h x Lx l) 175 x 338 x 196 mm
Weight: 3.9Kg
Classification: IIa

References

Monnal T20 KC030500
 Monnal T30 KC030600

Accessories:

Humidifier * KC030900
 Carrying bag KF007200
 Patient circuit KT600300
 Dust filter KY482500
 External battery * YR110301
 Software * KY666501

*Available mid-2009

Be prepared for
excellence